



Product Insert

PD Collagenase 100, PD Collagenase 800
Catalog # 011-3010/011-3020

Version Aug 2024

I. SUMMARY

Product Name	PD Collagenase 100 (PD 100), PD Collagenase 800 (PD 800)
Catalog Number	011-3010/011-3020
Grade	For Research Use Only
Stability (Expiry)	4 years from manufacture date
Storage (Lyophilized Cake)	≤ 2-8°C
Storage (Reconstituted Enzyme)	-20±5°C, 1 year without addition of other protease enzymes
Reconstitution Volume (mL)	Dependent on mass
Reconstitution Solutions	<ul style="list-style-type: none"> • RO/DI Water • Cell Culture Water • Water For Injection • HBSS Buffer • Lactated Ringers • Physiological Saline • RPMI • PBS (OK for immediate use of reconstituted enzyme, NOT for long term storage of reconstituted enzyme)
Reconstitution Time	Not less than 15 minutes
Animal Origin Statement	Bovine Free. Porcine gelatin peptone is used during the fermentation of <i>Clostridium histolyticum</i> but is largely removed by downstream processing. Non-mammalian gelatin peptone used in the final manufacturing of the final product.
Shipping	Enviro Ice Packs
Questions / Comments	<p>p. 317-917-3457</p> <p>e. feedback@vitacyte.com</p>

2. PRODUCT USE

2.1. Enzyme Reconstitution

While preparing for tissue digestion, equilibrate PD Collagenase to room temperature. PD Collagenase is supplied as a lyophilized powder. This powder may appear as a solid cake or clumps when first received. Vigorous shaking of the bottle or mechanical disruption with a laboratory spatula should quickly convert the material into a partially flowing powder. Weigh out the required amount of enzyme powder. The remaining enzyme may be resealed in the bottle and returned to storage at 2-8°C.

PD Collagenase may be reconstituted in a small volume of buffer or water and further diluted into the working buffer (suggest HBSS or a similar non-phosphate buffer) or added directly to the desired volume of working buffer. Allow the powder to rehydrate for a minimum of 15 minutes to ensure complete dissolution of the enzyme. Occasionally, invert the bottle to aid in the dissolution process. Enzyme denaturation may occur if the enzyme solution is vortexed or swirled excessively. The enzyme is lyophilized in a buffer containing calcium, so the initial reconstitution has sufficient calcium for enzyme stability. However, for optimal stability, the final working buffer for tissue dissociation should have at least 0.1 mM Ca²⁺ and contain no cation-chelating agents. The enzyme solution can be sterile filtered through 0.2 µm cellulose acetate or PES filter membranes without compromising enzyme potency. Surfactant-free cellulose acetate (SFCA) and PES filters from several major vendors were tested, and no measurable loss of collagenase's collagen degradation activity was observed. See details for preparing a mixture of the PD Collagenase 100 and 800 products in Appendix I.

2.2. Digestion Solution Preparation

See Appendices I & 2 for details on how to use the product.

3. PRODUCT DESCRIPTION

PD Collagenase 100 and 800 products are aseptically filled, lyophilized preparations of > 95% pure *Clostridium histolyticum* collagenases and purified neutral protease from *Paenibacillus polymyxa* with a composition detailed in the table below. These two products are equivalent to the DE Collagenase 100 and 800 products sold earlier by VitaCyte.

Components	1g PD Collagenase 100	1g PD Collagenase 800
<i>C. histolyticum</i> collagenase	55 mg	440 mg
BP Protease (Dispase™ equivalent enzyme)	18 mg	18 mg
Approximate mg polypeptide excipient	927 mg	542 mg

The minimally hygroscopic polypeptide excipient preserves enzyme stability during storage and adds convenience because users can weigh the precise amount of product needed immediately before use. These lyophilized products showed no loss of enzyme activity after stressing the lyophilized products by incubation at higher temperatures. These products have a five-year shelf life at 4°C.

The PD 100 product is formulated to have enzyme activities like those found in Worthington Type I collagenase, the first crude collagenase product used for cell isolation. The PD 800 product contains an eight-fold increase in collagenase mass. This formulation was defined by reverse engineering “good lots” of crude collagenase used for human islet isolation. The high collagenase-low protease activity mixture works well for many other cell types. This product contains the same amount of BP Protease found in the PD 100 product.

4. APPLICATION

- The PD Collagenase 100 and 800 alone or mixed in different ratios enable you to develop a continuum of collagenase to protease ratios for application-specific formulations.
- If the DE Collagenase product line was used in the past, see Appendix 1 to determine how to use the PD 100 and 800 products to prepare the enzyme formulations equivalent to DE 200, 400, or 600 products.
- Alternatively, these products can also be used to prepare purified enzyme mixtures to replace your existing product. To use this method, you will need to estimate the mass of collagenase used in the enzyme solution you prepare for cell isolation. See Appendix 2 for further details of this method.

5. STORAGE & STABILITY

PD Collagenase is stable for at least five years from the date of manufacture if stored as a lyophilized powder at $\leq 2-8^{\circ}\text{C}$. The product is shipped ambient but should be stored $\leq 2-8^{\circ}\text{C}$. Internal studies have shown the reconstituted enzyme is stable as a frozen solution at $-20\pm 5^{\circ}\text{C}$ for at least one year as long as no other protease enzymes had been added to the solution.

6. TROUBLESHOOTING

6.1. Many factors contribute to the successful isolation of cells from tissue and inadvertent oversight to any of these conditions may drastically reduce the yield and viability of target cell population. While far from a complete list, the guidance below may help identify commonly encountered problems. Contact VitaCyte if this guidance does not help resolve specific issues.

6.2. Prolonged or Incomplete Digestion may be caused by:

- Loss of enzyme potency (activity)
- Incomplete enzyme rehydration during reconstitution
- Inappropriate enzyme dilution
- Presence of enzyme inhibitors
- Low incubation temperature
- Inefficient digestion solution perfusion

6.3. Low Yield and/or Cell Viability

- Prolonged organ warm ischemia time
- Aggressive mechanical disruption
- Extended incubation time
- Elevated incubation temperature
- Inappropriate enzyme dilution


7. ADDITIONAL INFORMATION

7.1. Intended Use & Regulatory

PD Collagenases are for research use only.

7.2. Animal Origin

No bovine-derived animal products are used in any step of manufacturing PD Collagenases. Collagenase is purified from culture supernatants of *C. histolyticum* that contain porcine gelatin and pancreatic enzymes. Before lyophilization, a non-mammalian peptide excipient is added to the solution containing the purified enzyme mixture.

	Product Insert	
	PD Collagenase 100, PD Collagenase 800 <i>Catalog # 011-3010/011-3020</i>	Version Aug 2024

7.3. Manufacturing Summary

The purification processes use standard protein column chromatography, tangential flow filtration concentration, and diafiltration techniques. After characterization, the purified collagenases are sterile filtered in a qualified biosafety cabinet and aseptically dispensed by volume into amber bottles to contain 1 g dry weight of protein product. The final lyophilized product is then further characterized to confirm each batch meets established specification ranges.

7.4. Activity Assessment

Each lot of product is characterized for collagenase activity using the FALGPA peptide substrate¹ and neutral protease activity using succinyl casein substrate². The clostripain and trypsin-like activities are determined on the specific lot of collagenase used to prepare PD Collagenase products³.

7.5. Additional Considerations

In addition to the quality of the dissociation enzymes, additional factors impact the outcome of success of cell isolations including: the quality of the organ/tissue and experience of the cell isolation team. The team needs to assess many variables that affect islet recovery. These include but are not limited to the characteristics of the donor, transport of the organ/tissue, the cell isolation procedure, and subsequent cell culture.

7.6. Resources & Support

Further details on manufacturing, quality control testing and use of products are available at www.vitacyte.com or technical support at 317-917-3457.

7.7. References

- 1 Van Wart HE and Steinbrink DR. A continuous spectrophotometric assay for Clostridium histolyticum collagenase. Analytical Biochemistry 113 (1981); 356-65
- 2 Hatakeyma T, Kohzaki H, and Yamasaki N. A micro assay for proteases using succinylcasein as a substrate. Analytical Biochemistry 204 (1992); 181-184.
- 3 Mitchell WM and Harrington WF. Clostripain. Methods in Enzymology 19 (1970) 635-642
- 4 McCarthy RC, Breite AG, Green ML, Dwulet FE. Tissue dissociation enzymes for isolating human islets for transplantation: factors to consider in setting enzyme acceptance criteria. Transplantation 91 (2011) 137-45.
- 5 McCarthy RC, Green ML, Dwulet FE. Evolution of enzyme requirements for human islet isolation. OBM Transplantation 2 (2018) 024

8. APPENDICIES

-Appendix I

Refer to the table below to prepare products equivalent to the former DE Collagenase 200, 400, or 600 using PD Collagenase 100 & 800 products.



Product Insert

PD Collagenase 100, PD Collagenase 800
Catalog # 011-3010/011-3020

Version Aug 2024

Enzyme Formulation in reference to PD 100 product (Collagenase/Protease activity ratio)	Directions to prepare 1 mg equivalents of PD 200, PD 400, or PD 600 using PD 100 and PD-800	
	PD-Collagenase 100 (Cat # 011-1010)	PD Collagenase 800 (Cat # 011-1050)
DE 200 (2.6)	85.7% = 6 parts PD 100	14.3% = 1 part PD 800
DE 400 (5.2)	57.1% = 4 parts PD 100	42.9% = 3 parts PD 800
DE 600 (7.8)	28.6% = 2 parts PD 100	71.4% = 5 parts PD 800

The illustrative table shows how to prepare 100 mg of DE 200, DE 400, or DE 600

To check your calculations, add the mg of the PD 100 and PD 800. The answer should be the amount of DE 200, DE 400, or DE 600 products required for your cell isolation procedure.

Prepare the enzyme mixtures by weighing the appropriate mass of each product shortly before use—reconstitute and use within two hours. **DO NOT STORE POWDERED MIXTURES OF THE TWO PRODUCTS FOR RE-USE SINCE THE HOMOGENEITY OF THE NEW MIXTURE IS NOT KNOWN.**

	100 mg
DE 200	100 x 0.857 = 86 mg PD 100 + 100 x 0.143 = 14 mg PD 800
DE 400	100 x 0.571 = 57 mg PD 100 + 100 x 0.429 = 43 mg PD 800
DE 600	100 x 0.286 = 29 mg PD 100 + 100 x 0.714 = 71 mg PD 800

Appendix 2

The process described below applies the principles described in a [mechanistic model of enzyme-mediated cell isolation](#)^{4,5} to optimize collagenase and protease activities used in this process. The key elements of the model are summarized below.

- Cells are tethered to tissue by cell anchoring proteins within the extracellular matrix (ECM) which

are directly or indirectly associated with the collagen fibrils or fibers

- Collagen is the predominate protein in the ECM, providing protease resistant superstructure that protects/shields the ECM proteins from proteolytic degradation
- Proteases alone are unable to free cells from tissue, enzyme mixtures must contain a sufficient amount of collagenase's collagen degradation activity to degrade the collagen substrate
- Collagenase degrades collagen by the synergistic attack of class I (C1) and class II (C2) collagenase on individual tropocollagen molecules:
 - C1 is a processive protease that shaves of the collagen by moving along the molecule from the carboxy to the amino terminus
 - C2 is an endo protease that makes internal cuts on the collagen
- Collagenase disrupts the ECM, leading to exposure of protease sensitive sites on the ECM proteins
- Proteolytic cleavage of these exposed sites leads to proteases breaking the molecular tethers that hold cells to the ECM, leading to release of individual cells or cellular aggregates from tissue.


Four important principles are derived from this model:

- Purified collagenase has a narrow selectivity for collagen or gelatin (denatured collagen)
- Excess purified collagenase will have minimal adverse effect on cell viability or function because of its narrow to degrade collagen or gelatin
- The selection and dose of protease is the critical enzymatic activity that must be controlled to ensure the release of viable, functional cells from tissue
- The dose of neutral protease used for most cell isolations is excessive because it was determined by matching "good lots" of traditional collagenase products

The method described below applies these four principles to simplify the generation of a purified collagenase-protease enzyme mixture for your cell isolation procedure in 4 steps as described in the infographic below.


Step 1

Estimate collagenase mass in product



Step 2

Calculate mass of PD Collagenase needed to prepare an equivalent collagenase mass $A = B$ ($B = \text{PD Collagenase}$)




Step 3

Prepare 3 enzyme solution (1, 2, 3) each containing equivalent amounts of collagenase (B), but different amounts of protease activity: +++ to ++ to +



Step 4

Review yields of functional cells, selected best PD Collagenase mixture, Reconfirm results in additional isolations



Step 1: Find the collagenase mass in the product you plan or currently use for isolating cells. Roche Liberase research products are sold in 5 or 50 mg bottles of collagenase. To determine the mg of collagenase in Worthington Types 1, 2, 3, or 4 collagenase products, divide the dry weight of collagenase to prepare your enzyme solution by 16.67 (these product contain about 6% collagenase). For other products, contact VitaCyte technical support.

For reference, the table below shows the mg of collagenase in the enzyme solutions used to isolate the specific cells below. A reference is provided to show the method used for cell isolation. These values were derived from consulting the life science protocol websites listed under Applications > Other Cells tab on the website or by performing literature searches.



Product Insert

PD Collagenase 100, PD Collagenase 800
Catalog # 011-3010/011-3020

Version Aug 2024

Tissue	Species	Cell isolated	mg/ml collagenase enzyme solution	Reference
Adipose parametrial fat	Rat	ADSVF	3.0	https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4684382/
Heart	Adult rat or mouse	Cardiomyocytes	0.06 for mouse 1.6 mg for rat	https://pubmed.ncbi.nlm.nih.gov/35813022/
Venous	Adult humans umbilical tissue	HUVECs	0.07	https://pubmed.ncbi.nlm.nih.gov/18978951/
Kidney	Adult mice	Renal tubular epithelial cells	0.25	https://pubmed.ncbi.nlm.nih.gov/29985358/
Skeletal muscle	Rats	Satellite cells	0.15	https://pubmed.ncbi.nlm.nih.gov/23542587/

Step 2: Determine the amount of collagenase required for your isolation procedure. This amount is the mass you use to prepare the enzyme solution in the cell isolation process. Make the appropriate calculations to determine the mass of collagenase in the enzyme solution you use to isolate cells.

Step 3: Prepare three enzyme solutions for lot testing using the PD Collagenase 100, PD 800 products and a 50:50 mixture of PD 100:PD 800. The amount of collagenase, protease and excipient for these 3 products are listed in the table below.

Product Components	mg of each component in 1 g product		
	PD Collagenase 100	50:50 mixture of PD 100:PD 800	PD Collagenase 800
≈ mg <i>C. histolyticum</i> collagenase	55	248	440
mg BP Protease (Dispase™ equivalent enzyme)	18	18	18
≈ mg of stabilizing polypeptide excipient	927	734	542
Percent collagenase per dry weight product	5.5	24.8	44.0

For illustrative purposes, assume you want to make an enzyme solution with 10 mg of collagenase (note, collagenase mass ≠ dry weight of traditional collagenase product) in 100 mL. You can assume the collagenase is in excess if you successfully isolate your cells of interest most of the time.

Three enzyme solutions are prepared, each containing 10 mg of collagenase. The amount is determined by dividing 10 by the percentage of collagenase in each product.

- For PD 100: $10 \text{ mg} / 0.055 = 181.8 \text{ mg}$
- For 50:50 mixture of PD100:PD800: $10 \text{ mg} / 0.248 = 40.3 \text{ mg}$
- For PD 800: $10 \text{ mg} / 0.44 = 22.7 \text{ mg}$

The table below shows that each enzyme solution contains the same amount of collagenase (10 mg in 100 mL or 100 ug/mL of enzyme solution) but differing amounts of neutral protease.

Enzymes	PD Collagenase 100	50:50 mixture of PD 100:PD 800	PD Collagenase 800
Total mg collagenase in enzyme solution	10	10	10
Total mg BP Protease in the enzyme solution	3.27	0.72	0.41
ug BP Protease per mL enzyme solution	32.7	7.2	4.1
Neutral protease activity (NPA U/mL)	4578	1008	574

Step 4: Determine if the VitaCyte's products provides comparable results. You may find slower digestion times as the neutral protease activity is decreased. This is expected from the hypothetical model of enzyme mediated tissue digestion. You will need to determine if extending the digestion time provides other benefits to recovered cell population.

If the results are comparable to those obtained with your current lot of product, congratulations! You now know more about the enzyme formulation required for isolating your cell of interest.

If you do not obtain equivalent results with any of the enzyme mixtures above when compared to your current lot of collagenase, the enzyme formulation may need to adjusted or supplemented with another protease, Clostripain, as described in the package insert. Clostripain may be needed because it has a complementary enzyme activity when compared to BP Protease. Clostripain is a trypsin-like enzyme that cuts at different regions of the protein than BP Protease.

The effort required to assess the performance of the three PD Collagenase enzyme solutions above will not change the effort you presently take to qualify new lots of traditional crude or enriched collagenase products. The benefits of adopting the PD Collagenase enzyme in place of traditional collagenase are summarized in the table below.



Product Insert

PD Collagenase 100, PD Collagenase 800
Catalog # 011-3010/011-3020

Version Aug 2024

	Lot Qualification	PD Optimization
Effort required	Equivalent	Equivalent
Knowledge gained	None, no knowledge of enzyme composition	Essential: enzyme composition defined
Lot pre-qualification	No change, must prequalify new lots	Once defined, no need to prequalify future lots
Ability to modify formulation	None	Yes
Shelf stability	?	4 years

If you have any questions on using enzymes for cell isolation or recovery, contact VitaCyte Technical Support.